



Bruker BioSpin MRI

ParaVision 360 V3.6

# Shim Calculation

- Protocol "T2star\_map\_MGE"

# Table of Contents

1. Basic Shim Results .....	2
2. Acquisition Information .....	3
3. System Information .....	4

## 1. Basic Shim Results

Status of shim calculation

Succeeded.

**Table 1. Used Pixels**



Used Field Values	4,790
Effective Volume	285.506ul

**Table 2. Basic Shim Vector**

Coil	Map Shim /%	New Shim /%	New Shim /A	Term	Unit
Z0	0.000	0.093	0.001	-582.4	Hz (Z0)
Z	-10.499	-5.380	-0.092	-123.7	Hz / cm**1 (Z)
Z2	0.000	5.154	0.258	184.3	Hz / cm**2 (Z2)
Z3	0.000	-0.324	-0.016	-1.3	Hz / cm**1 (Z)
X	2.782	2.753	0.047	67.5	Hz / cm**1 (X)
Y	-5.198	-7.591	-0.130	-184.0	Hz / cm**1 (Y)
ZX	0.000	0.338	0.017	13.0	Hz / cm**2 (XZ)
ZY	0.000	0.119	0.006	4.7	Hz / cm**2 (YZ)
2XY	0.000	-0.890	-0.044	-22.8	Hz / cm**2 (XY)
X2-Y2	0.000	1.139	0.057	27.3	Hz / cm**2 ((X2-Y2))

**Table 3. Basic Shim Statistics**

Property	Map	Basic Shim (est.)
Mean/Hz	-8.5	0.0
Standard Deviation /Hz	16.6	8.8
Min/Hz	-103.5	-107.9
Max/Hz	84.7	98.0
Absolute Dev./Hz	12.1	4.9
Outliers Low/High	37 / 12	79 / 72

**Figure 1. Histogramm of measured Field Values****Figure 2. Histogramm of Predicted Field Values**

## 2. Acquisition Information

Study Directory:

20240725\_090212\_std\_PV360\_3\_6\_94T\_protocols\_1\_1

**Table 4. Field Map**

Property	Value
Expno	19
Procno	1
Scan Cradle position	0
Map Cradle position	0.002

**Table 5. Basic Shim Area**

Property	Value
Shape	Ellipsoid_In_GobjShape
Extent +1; 1st dir	10.175 mm
Extent +2; 2nd dir	5.798 mm
Extent +3; 3rd dir	9.451 mm
Margin +1; 1st dir	0 mm
Margin +2; 2nd dir	0 mm
Margin +3; 3rd dir	0 mm

### 3. System Information

**Table 6. System Information**

Electronic System	AVANCE NEO
Field Strength	9.404 T
Frequency	400.395 MHz

**Table 7. Gradient Hardware Information**

Gradient Coil	B-GA105S HP GRAD SYSTEM BC94/17 MAXWELL (BFR_W149640_0021)
Gradient Amplifier Type	BGA
Gradient Power Supplies	IECO_XPA

Shim chain identifier:

B-GA105S HP GRAD SYSTEM BC94/17 MAXWELL (BFR W149640/0021).222714468

Parameter that determine the shim chain identifier are marked with a (\*).

**Table 8. Shim Hardware Information**

Shim Coil (*)	B-GA105S HP GRAD SYSTEM BC94/17 MAXWELL (BFR_W149640_0021)
Shim Power Supply	BSPS_8155
Total SPS ouptut channels (*)	8
Shim Unit Type (*)	GTCUBE
Max Shim Value / Hardware units (*)	131,070

**Table 9. Shim Power Supply Output Channels**

SPS Chan To Coil List (*)	SPS Chan To Coil I <sub>max</sub> (*)	SPS Chan To Coil U <sub>max</sub> (*)
Shim_Z2	5 A	16 V
Not Connected	5 A	16 V
Shim_YZ	5 A	16 V
Shim_XZ	5 A	16 V
Shim_X2_Y2	5 A	16 V
Shim_2XY	5 A	16 V
Shim_Z3	5 A	16 V
Not Connected	5 A	16 V

**Table 10. Parameter that determine the shim chain identifier**

Shim Checksum Component Names (*)	Shim Checksum Component Values (*)
GTCUBE	ShimUnit=7
shimMaxVal=131070	shimMaxVal=131070
nSPS=8	nSPS=8
Shim_Z2	SPS_ch2sh[0]=2
SPS_vmax[0]=16	SPS_vmax[0]=16
SPS_imax[0]=5000	SPS_imax[0]=5000
Shim_Invalid	SPS_ch2sh[1]=64
SPS_vmax[1]=16	SPS_vmax[1]=16
SPS_imax[1]=5000	SPS_imax[1]=5000
Shim_YZ	SPS_ch2sh[2]=8
SPS_vmax[2]=16	SPS_vmax[2]=16
SPS_imax[2]=5000	SPS_imax[2]=5000
Shim_XZ	SPS_ch2sh[3]=7
SPS_vmax[3]=16	SPS_vmax[3]=16
SPS_imax[3]=5000	SPS_imax[3]=5000
Shim_X2_Y2	SPS_ch2sh[4]=10
SPS_vmax[4]=16	SPS_vmax[4]=16
SPS_imax[4]=5000	SPS_imax[4]=5000
Shim_2XY	SPS_ch2sh[5]=9
SPS_vmax[5]=16	SPS_vmax[5]=16
SPS_imax[5]=5000	SPS_imax[5]=5000
Shim_Z3	SPS_ch2sh[6]=3
SPS_vmax[6]=16	SPS_vmax[6]=16
SPS_imax[6]=5000	SPS_imax[6]=5000
Shim_Invalid	SPS_ch2sh[7]=64
SPS_vmax[7]=16	SPS_vmax[7]=16
SPS_imax[7]=5000	SPS_imax[7]=5000
shimStatus=BFR_W149640_0021	shimStatus=BFR_W149640_0021
nGPS=3	nGPS=3

Shim Checksum Component Names (*)	Shim Checksum Component Values (*)
IECO_XPA	GPS_type[1]=15
GPS_imax[1]=150000	GPS_imax[1]=150000
IECO_XPA	GPS_type[2]=15
GPS_imax[2]=150000	GPS_imax[2]=150000
IECO_XPA	GPS_type[4]=15
GPS_imax[4]=150000	GPS_imax[4]=150000
gradientStatus=BFR_W149640_0021	gradientStatus=BFR_W149640_0021

**Table 11. Shim Coil Properties**

Coil Name	Coil Identifier	SPS Index	SPS I <sub>max</sub> /A
Z0	Ch1: Spf17 FRED	0	1
Z	Ch2: Spf3 Gradient Z	0	1.705
Z2	Ch3: Spf4 SPS1,1 CoilPins=B1-B2 CoilId=1	1	5
Z3	Ch4: Spf9 SPS1,7 CoilPins=B4-B5 CoilId=6	1	5
X	Ch5: Spf1 Gradient X	0	1.713
Y	Ch6: Spf2 Gradient Y	0	1.713
ZX	Ch7: Spf5 SPS1,4 CoilPins=C4-C5 CoilId=2	1	5
ZY	Ch8: Spf6 SPS1,3 CoilPins=D1-D2 CoilId=3	1	5
2XY	Ch9: Spf7 SPS1,6 CoilPins=F1-F2 CoilId=4	1	5
X2-Y2	Ch10: Spf8 SPS1,5 CoilPins=E4-E5 CoilId=5	1	5

**Table 12. Gradient Hardware Components**

Gradient Power Supply (*)	IECO_XPA
Max Current (*)	150 A